PERIPHERALLY INSERTED CENTRAL
CATHETER (PICC)

INSTRUCTIONS FOR USE:

DEVICE DESCRIPTION

- The PICC product family is made from biocompatible radiopaque polyurethane material. The PICC is offered in different diameters, and number of lumens. They are designed to facilitate frequent blood sampling, delivery of medications, nutrition, blood products, central venous pressure monitoring and injection of contrast media solutions.

DEVICE MATERIAL COMPOSITION

- The materials and their composition that make up the PICC catheter are: PVC (11-28%, 1.38 g – 3.87 g), Acetal (15-18%, 1.23 g – 2.46 g), ABS (10-11%, 0.92 g – 1.46 g), Stainless Steel (5-9%, 0.71 g), and Polyurethane (26-51%, 2.10 g – 6.38 g). The portion of the catheter that remains in the vessel is made of 100% Polyurethane.

INTENDED PURPOSE

- The PICC catheter is a medical device intended for adult patients that need short or long term (less than or greater than 30 days) peripheral access to the central venous system.

INDICATIONS FOR USE

The PICC catheter is indicated for adult critically ill patients where short or long term peripheral access to the central venous system is required for:

- Intravenous therapy
- Infusion of vesicants, irritants, antibiotics, and medications
- Administration of nutrients for partial or total parenteral nutrition
- Frequent blood sampling collection
- Injection of contrast media for diagnostic imaging exams
- Central venous pressure monitoring

IMPORTANT INFORMATION FOR POWER
INJECTION

- Use only lumens marked “Power Injectable” for power injection of contrast media.

CAUTION: Use of lumens not marked “Power Injectable” for power injection of contrast media may cause failure of the catheter.

- Contrast media should be warmed to body temperature prior to power injection.

WARNING: Failure to warm contrast to body temperature prior to power injection may result in catheter failure.

- Vigorously flush the PICC catheter using a 10cc or larger syringe and sterile normal saline prior to and immediately following the completion of power injection studies. This will ensure the patency of the catheter and prevent damage to the catheter. Resistance to flushing may indicate partial or complete catheter occlusion. Do not proceed with power injection study until occlusion has been cleared.

WARNING: Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.

WARNING: Do not exceed the maximum indicated flow rate printed on the catheter.

WARNING: Power injector machine pressure limiting feature may not prevent over pressurization of an occluded catheter.

WARNING: Exceeding the maximum indicated flow rate may result in catheter failure and/or catheter tip displacement.

WARNING: The indication of power injection of contrast media implies the catheter's ability to withstand the procedure but does not imply appropriateness of the procedure for a particular patient. A suitably trained clinician is responsible for evaluating the health status of a patient as it pertains to a power injection procedure.

CONTRAINDICATIONS

- The presence of infection, bacteremia, or septicemia is known or suspected.
- The patient's anatomy is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- There has been past irradiation of prospective insertion site.
- There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- There are local tissue factors that may prevent proper device stabilization and/or access. Or skin conditions that prevent insertion through an unaffected area.

POSSIBLE COMPLICATIONS

- Accidental arterial access
- Air embolism
- Bleeding
- Brachial plexus injury
- Cardiac arrhythmia
- Cardiac tamponade

- Catheter breakage
- Catheter embolism
- Catheter erosion through the skin
- Catheter occlusion
- Catheter related blood stream infection (CRBSI)
- Catheter related sepsis
- Central line associated blood stream infection (CLABSI)
- Endocarditis
- Exit site infection
- Exit site necrosis
- Extravasation
- Fibrin sheath formation
- Hematoma
- Intolerance reaction to implanted device
- Myocardial erosion
- Phlebitis
- Spontaneous catheter tip malposition or retraction
- Thromboembolism
- Venous thrombosis
- Ventricular thrombosis
- Vessel injury
- Viscus injury
- Risks normally associated with local or general anesthesia, surgery, and post-operative recovery

WARNINGS

- Insertion of these devices requires training consistent with the Standards of Practice. Before attempting the insertion, ensure that you are familiar with the complications and their emergency treatment should any of them occur.
- In the rare event that a hub or connector separates from any component during insertion or use, take all necessary steps and precautions to prevent blood loss or air embolism and remove the catheter.
- Do not advance the guidewire or catheter if unusual resistance is encountered.
- Do not insert or withdraw the guidewire quickly or forcibly from any component. The wire may break or unravel. If the guidewire becomes damaged, the introducer needle or sheath/dilator and guidewire must be removed together.
- Federal Law (USA) restricts this device to sale by or on the order of a physician.
- This catheter is for single use only, do not reuse. Reuse of catheter may result in cross-contamination, plastic degradation, catheter rupture, etc.
- Do not re-sterilize the catheter or accessories by any method.
- The manufacturer shall not be liable for any damages caused by reuse or re-sterilization of this catheter or accessories provided.

- Contents sterile and non-pyrogenic in unopened, undamaged package.
- Do not use catheter or accessories if package is opened or damaged.
- Do not use catheter or accessories if any sign of product damage is visible.
- Do not perform injection of contrast media in a patient receiving hemodialysis treatment.
- Do not use in pediatric patients.
- Do not overtighten luer access devices
- Manufacturer does not warrant the use of any accessories that are not provided with catheter.
- The stylet used in the positioning of the PICC is composed of stainless steel and contains < 2 % cobalt. Cobalt is considered a carcinogen under (EC) No. 1907/2006. However, assessment of exposure found that the increased cancer risk associated with the stylet's use is insignificant as defined under ISO 10993-17.

CATHETER PRECAUTIONS

- Small syringes will generate excessive pressure and may damage the catheter. Ten (10)cc or larger syringes are recommended.
- Do not use sharp instruments near the extension lines or catheter lumen.
- Do not use scissors to remove dressing.
- Catheter will be damaged if clamps other than what is provided with this kit are used.
- Clamping of the tubing repeatedly in the same location will weaken tubing. Avoid clamping near the luer(s) and hub of the catheter.
- Examine catheter lumen and extension(s) before and after each infusion for damage.
- To prevent accidents, assure the security of all caps and connections prior to and between uses.
- Use only Luer Lock (threaded) Connectors with this catheter.
- Repeated over tightening of luer lock connections, syringes, and caps will reduce connector life and could lead to potential connector failure.
- Confirm catheter tip position by x-ray, fluoroscopy and/or catheter tip confirmation technology prior to use. Monitor tip placement routinely per institution policy.
- The maximum pressure safety cutoff setting of the power injector utilized should not exceed the pressure of the catheter used as shown in the table below.


- Power injector machine pressure limiting feature may not prevent over pressurization of an occluded catheter.

French/ Lumen(s)	Rated occluded catheter burst pressure (Psi)
3Fr/ Single	352 Psi
4Fr/ Single	341 Psi
4Fr/ Dual	327 Psi
5Fr/ Single	330 Psi
5Fr/ Dual	308 Psi
6Fr/ Dual	291 Psi
6Fr/ Triple	301 Psi

INSERTION SITES

- The basilic, brachial, or cephalic vein may be catheterized. The basilic vein above antecubital fossa is the preferred site.
- CAUTION:** For PICC catheter with a reverse-taper design. Placement of a larger catheter at or below the antecubital fossa may result in an increased incidence of phlebitis. Placement of the PICC catheter above the antecubital fossa is recommended.

Basilic Vein

DIRECTIONS FOR MODIFIED SELDINGER
INSERTION

- Read instructions carefully before using this device. The catheter should be inserted, manipulated, and removed by a qualified, licensed physician or other qualified health care professional under the direction of a physician.
 - The medical techniques and procedures described in these instructions for use do not represent all medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.
 - Use standard hospital protocols when applicable.
- PRIOR TO PLACEMENT**
- Identify Insertion Site and Vein, considering the following variables:
- Patient diagnosis

- Age and size of patient
- Unusual anatomical variables
- Type and purpose of IV therapy
- Anticipated dwell time of catheter

1. Apply tourniquet to arm above anticipated insertion site.
2. Select vein based on assessment.
3. Release tourniquet.

PREPARE CATHETER

- Preflush catheter, sideport adapter, and needleless access port(s).
- Attach saline filled syringe to luer of sideport adapter and flush adapter and catheter. Clamp sideport extension and remove syringe. If using double lumen catheter, attach needleless access port to remaining extension. Attach saline filled syringe to the needleless access port and completely flush catheter lumen. Remove syringe from needleless access port prior to clamping extension.

CAUTION : Never close clamp on catheter stylet; stylet and catheter damage may result.

CAUTION : The needleless access port should not be used with needles, blunt cannula, or other non-luer connectors, or luer connectors with visible defects. If needle access is attempted, the needleless access port must be replaced immediately. Do not exceed 100 actuations.

INSERTION

- Strict sterile technique must be used during insertion, maintenance, and catheter removal procedures. Provide a sterile operative field. Use sterile drapes, instruments, and accessories. Perform surgical scrub. Wear gown, cap, sterile gloves, and mask.
 - Apply tourniquet to arm above anticipated insertion site to distend the vein.
 - Insert the introducer needle into the target vein. Attach a syringe and aspirate to ensure blood return. Release tourniquet.
 - Removing the syringe, place thumb over the end of the needle to prevent blood loss or air embolism. Draw the flexible end of marked .018" guidewire back into advancer so that only the end of the guidewire is visible. Insert the advancer's distal end into the needle hub. Advance guidewire with forward motion into and past the needle hub into the target vein.
- CAUTION :** The length of the wire inserted is determined by the size of the patient. Monitor patient for arrhythmia throughout this procedure. The patient should be placed on a cardiac monitor during this procedure. Cardiac arrhythmias may result if guidewire is

allowed to pass into the right atrium. The guidewire should be held securely during this procedure.

Note: For alternate insertion method, see Directions for Seldinger Insertion Section.

- Remove needle, leaving guidewire in the target vein. Thread sheath/dilator over the proximal end of the guidewire into target vein. Advance the guidewire until it reaches the cavoatrial junction. Once the guidewire is in place, measure the depth of the guidewire by reading the markings on the wire. Remove the guidewire leaving the sheath and dilator in the vein.

CAUTION: DO NOT bend the sheath/dilator during insertion as bending will cause the sheath to prematurely tear. Hold sheath/dilator close to the tip (approximately 3cm from tip) when initially inserting through the skin surface. To progress the sheath/dilator towards the vein, regrasp the sheath/dilator a few centimeters (approximately 5cm) above the original grasp location and push down on the sheath/dilator. Repeat procedure until sheath/dilator is fully inserted.

CAUTION: Never leave sheath in place as an indwelling catheter. Damage to the vein will occur.

- Loosen locking collar of sideport and withdraw stylet back beyond the point where the catheter is to be trimmed by at least ¼ inch (1 cm), Cut catheter to length determined by marked guidewire.

CAUTION: Never attempt to cut stylet.

CAUTION: Always withdraw stylet back beyond the tip of the catheter prior to insertion.

- Once proper catheter length and stylet position has been achieved, tighten locking collar to keep stylet in place.
- Remove dilator from sheath.
- Insert distal tip of catheter into and through the sheath until catheter tip is correctly positioned in the target vein.
- Remove the tear-away sheath by slowly pulling it out of the vessel while simultaneously splitting the sheath by grasping the tabs and pulling them apart (a slight twisting motion may be helpful)

CAUTION: Do not pull apart the portion of the sheath that remains in the vessel. To avoid vessel damage, pull back the sheath as far as possible and tear the sheath only a few centimeters at a time.

- Make any adjustments to catheter under fluoroscopy. The distal tip should be positioned at the level of the cavoatrial junction.
 - Remove the stylet by applying gentle pressure with one hand above the insertion site while grasping the stylet with the other hand and slowly pulling back with a constant motion.
- CAUTION:** Do not use excessive force to remove stylet. Slightly loosen stylet collar to reduce resistance during removal.
- Remove sideport adapter and replace with needleless access port. Attach saline filled syringe to needleless access port, aspirate lumen and then irrigate with saline. Remove syringe prior to clamping extension.
- CAUTION:** If difficulty and/or bunching of the catheter lumen are experienced while removing the stylet, additional flushing of the catheter may be helpful. The catheter may need to be repositioned to allow for removal of the stylet.
- CAUTION:** Do not attempt to reinsert stylet once it has been withdrawn.
- CAUTION:** Never leave stylet in place after catheter insertion; injury may occur. Remove both stylet and sideport adapter after insertion.
- Attach syringe(s) to extension(s) and open clamp(s). Blood should aspirate easily. If excessive resistance to blood aspiration is experienced, the catheter may need to be repositioned to obtain adequate flow.
 - Once adequate aspiration has been achieved, lumen(s) should be irrigated with saline filled syringe(s). Clamp(s) should be open for this procedure.
- CAUTION:** Small syringes will generate excessive pressure and may damage the catheter. Ten (10)cc or larger syringes are recommended.
- Remove the syringe(s) and close extension clamp(s). Avoid air embolism by keeping catheter tubing clamped at all times when not in use and by aspirating then irrigating the catheter with saline prior to each use. With each change in tubing connections, purge air from the catheter and all connecting tubing and caps.
- CAUTION:** Do not clamp the lumen portion of the catheter. Clamp only the extension(s). Do not use the serrated forceps, use only the in-line clamp(s) provided.
- Confirm and document proper tip placement with fluoroscopy prior to use. The

distal tip should be positioned at the level of the cavoatrial junction.

CAUTION: Failure to verify catheter placement may result in serious trauma or fatal complications.

Note: If there is no blood return, verify catheter position before use.

CATHETER SECUREMENT AND WOUND DRESSING:

- The insertion site and external portion of the catheter should always be covered with a protective dressing.
- Cover the exit site with an occlusive dressing according to the facility policy.
- Record catheter length, catheter lot number, and tip position on patient's chart.

DIRECTIONS FOR SELDINGER INSERTION
- Follow directions for Modified Seldinger Insertion, up to INSERTION step (bullet) #5. - Remove needle, leaving guidewire in the targeted vein. Advance the guidewire until it reaches the cavoatrial junction. Once the guidewire is in place, measure the depth of the guidewire by reading the markings on the wire. - Remove stylet from catheter. Failure to do so will result in difficulty advancing catheter over the wire - Cut catheter to length determined by marked guidewire. - Insert proximal end of wire into distal tip of catheter lumen. Feed catheter lumen into the vessel following the guidewire. Advance catheter lumen along the guidewire until the distal tip is correctly positioned in the target vein. The distal tip should be positioned at the level of the cavoatrial junction.
CAUTION: A skin nick may be required to feed the catheter smoothly into the vessel. - Make any adjustments to catheter under fluoroscopy. The distal tip should be positioned at the level of the cavoatrial junction. CAUTION: Do not clamp the lumen portion of the catheter. Clamp only the extension(s). Do not use serrated forceps; use only in-line clamps provided. - Remove the wire from the catheter. Remove by applying gentle pressure with one hand above the insertion site while grasping the 130cm wire with the other hand and pulling slowly back with a constant motion. - Follow directions for Modified Seldinger Insertion, from INSERTION step (bullet) # 14 on.

POWER INJECTION PROCEDURE
Remove the injection/needleless cap from the PICC catheter.
- Using a 10cc or larger syringe(s), aspirate catheter lumen(s) to assure patency and remove heparin. Discard syringe(s). - Attach a 10cc or larger syringe filled with sterile normal saline and vigorously flush the catheter with the full 10cc of sterile normal saline.
WARNING: Failure to ensure patency of the catheter prior to power injection studies may result in catheter failure.
- Detach syringe. - Attach the power injection device to the PICC catheter per manufacturer's recommendations.
WARNING: Always use connector tubing between power injector syringe and catheter. Do not attempt to connect power injector syringe directly to the catheter. Damage may result.
Complete power injection study taking care not to exceed the flow rate limits.
WARNING: Exceeding the maximum indicated flow rate may result in catheter failure and/or catheter tip displacement.
- Disconnect the power injection device. - Flush the PICC catheter with 10cc of sterile normal saline, using a 10cc or larger syringe. If double lumen catheter, flush both lumens. - Replace the injection/needleless cap on the PICC Catheter.

French/ Lumen(s)	Rated Power Injection Flow Rate (Psi)
3Fr/ Single	0.5 mL/sec
4Fr/ Single	5 mL/sec
4Fr/ Dual	2 mL/sec
5Fr/ Single	5 mL/sec
5Fr/ Dual	
6Fr/ Dual	
6Fr/ Triple	

INFUSION
- Before infusion begins all connections should be examined carefully. - Frequent visual inspection should be conducted to detect leaks to prevent blood loss or air embolism. - If a leak is found, the catheter should be clamped immediately and replaced. - Necessary remedial action must be taken prior to the continuation of the treatment.
Note: Excessive blood loss may lead to patient shock.

CATHETER MAINTENANCE
Dressing Changes-A dressing should always cover the insertion site . The dressing should be changed per institutional policy or any

time the dressing becomes soiled, wet, or non-occlusive.

Note: When using alcohol or alcohol containing antiseptics with the PICC, care should be taken to avoid prolonged or excessive contact. Solutions should be allowed to completely dry before applying an occlusive dressing. Chlorhexidine gluconate and/or povidone iodine are the suggested antiseptics to use.

WARNING : Alcohol should not be used to soak or declot the PICC because alcohol is known to degrade polyurethane catheters over time with repeated and prolonged exposure.

WARNING : Acetone and polyethylene glycol containing ointments should not be used with the PICC, as these may cause failure of the device.

Note: During all dressing changes the external length of the catheter should be measured to determine if catheter migration has occurred. Periodically confirm catheter placement and tip location by imaging method.

- Flushing and Heparinization – Follow institutional policy for flushing frequency and heparin concentration.
- The catheter should be flushed with normal saline prior to drug administration to remove heparin solution.
- After drug administration each lumen should be flushed again with normal saline and then locked with heparin to maintain patency.

Injection Caps - Injection cap(s) or needleless access port(s) should be changed per institutional policy.

CATHETER PERFORMANCE
Occluded/Partially Occluded Catheter- If resistance is encountered to aspirating or flushing, the lumen may be partially or completely occluded.
WARNING : Do not flush against resistance.
If the lumen will neither aspirate nor flush, and it has been determined that the catheter is occluded with blood, follow institutional declothing procedure.
Infection
CAUTION : Due to risk of exposure to HIV or other blood borne pathogens, health care professionals should always use Universal Blood and Body Fluid Precautions in the care of all patients.
Sterile technique should always be strictly adhered to.

EXPECTED LIFETIME
- The PICC is a medical device indicated for short or long term, less than or equal to 12-months, however, the lifetime of the PICC is affected by many factors related to the use of the device, such as: - Incorrect placement and insertion sites other than indicated in the instructions for use. - Prolong or excessive contact with antiseptic containing alcohol or other alcohol solutions not indicated in the instructions for use which may cause damage or shorten the life of the catheter. - Using small syringes less than 10mL generate excessive pressure and damage the catheter. - Clamping of the tubing repeatedly in the same location will weaken the extension tubing. - Using sharp instruments near the extension lines or catheter lumen. - Using scissors when removing or changing wound dressings. - Overtightening of luer access devices

CATHETER REMOVAL
WARNING : Only a clinician familiar with the appropriate techniques should attempt the following procedures.
CAUTION: Always review facility protocol, potential complications and their treatment, warnings, and precautions prior to catheter removal.
- Wash hands gather equipment. - Remove old dressing and inspect insertion site for redness, tenderness, and drainage. - Grasp catheter near insertion site and using a slow steady motion, remove catheter from vein. - If resistance is felt - STOP. Retape the catheter and apply a warm compress to the extremity for 20-30 minutes. - Resume removal procedure. If further difficulty is encountered, follow institutional policy for further intervention. - Apply pressure, if necessary, until bleeding stops and dress site following institutional policy.
Note: Inspect catheter and measure length. It must be equal to baseline measurement taken when the catheter was inserted.

SAFETY DISPOSAL
Discard all procedural components according to hospital protocol in an appropriate

biohazard and/or sharp disposal container following applicable regulations.

ADDITIONAL INFORMATION
- The catheter should be inserted, manipulated, and removed only by a qualified, licensed physician or other qualified health care professional under the direction of a physician. - The PICC catheter is intended only for adult patients. - The Summary of Safety and Clinical Performance (SSCP) is available in the European database on medical devices (EUDAMED) linked to the Basic UDI-DI: 087660400PICC01001KX For the SSCP visit: www.ec.europa.eu/tools/eudamed - For e-IFU visit: www.hlic.net/products/picc-lines - For any serious incident that has occurred in relation to the device, contact: EU Auth. Rep.: emergoeurope@ul.com Manufacturer: customerservice@hlic.net
DEVICES IN CONNECTION WITH THE PICC DEVICE

The PICC catheter interacts during a standard procedure with the following devices:

- The stylet during the catheter insertion to control the movement of the catheter as it is inserted through the vein to confirm safe passage.
- The peelable tearaway introducer during the catheter insertion guides the catheter through the vein until the catheter tip is correctly positioned in the target vein.
- The grip-lok securement after catheter placement stabilized the catheter in a fixed position for protection and safety.
- The syringe after catheter insertion flushed the catheter and extensions with saline fluid.
- The needleless connector after catheter insertion connects the catheter to other external devices.

CLINICAL BENEFITS
- Preserve vessel health from vesicant and irritant medication. - Reduce the need for repeated venous catheter insertion. - Provide venous access for IV nutrition and IV therapies. - Provide venous access for central venous pressure monitoring. - Allow patients with chronic medical conditions to receive long term treatment - Allow patients to receive administration of compatible and non-compatible medications (precipitates). - Allow physicians to provide rapid determination on applicable fluid therapy to patients.

• Allow CT diagnostic imaging of IV contrast media

WARRANTY
Health Line International Corp. WARRANTS THAT THIS PRODUCT WAS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS. PATIENT CONDITION, CLINICAL TREATMENT, AND PRODUCT MAINTENANCE MAY AFFECT THE PERFORMANCE OF THIS PRODUCT. USE OF THIS PRODUCT SHOULD BE IN ACCORDANCE WITH THE INSTRUCTIONS PROVIDED AND AS DIRECTED BY THE PRESCRIBING PHYSICIAN.
FR
 HEALTH LINE MEDICAL PRODUCTS
CATHÉTER CENTRAL À INSERTION PÉRIPHÉRIQUE (PICC) MODE D’EMPLOI :
DESCRIPTION DU DISPOSITIF

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Health Line International Corporation