# **ProV Midline ST**


Midline Catheters with V-Precision Tip





**ProV Midline ST** is a Midline catheter intended for short term (less than 30 days) peripheral access to the venous system for the purpose of intravenous therapy, medicines, and blood products. The smooth, atraumatic tipped cannula allows for smooth insertion and helps to avoid trauma to the vessel.

A bridge between a short PIVC and a traditional Midline, Health Line's **ProV Midline ST** is a short midline catheter, a suitable choice for patients requiring longer cannulation than achieved by a short PIVC, and whose treatment regiment includes frequent blood draws, power injection cycles, and multiple infusions within 30 days.



## **Features and Benefits**

Power Injectable: Allows for the power injection of contrast media.

**'V' Tipped Catheter Design:** Aides in the transition from wire to vessels while helping promote consistent asperation and flow rates.

Seldinger Technique: Facilitates a simple and straightforward, OTW insertion.

**Multiple Sizes**: The ProV Midline ST is available in multiple sizes to accommodate various vessel sizes and aid in patient comfort.

Multiple Uses: Designed to be used for infusion, intravenous therapy, and blood sampling.

### **Health Line Prov Midline ST Catheters**

#### 3FR Kit Range

Catalogue No.	Description	Box Qty	Carton Qty
A130S-03106	3Fr x 06cm Single Lumen Midline Seldinger Kit	5	20
A130S-03108	3Fr x 08cm Single Lumen Midline Seldinger Kit	5	20
A130S-03110	3Fr x 10cm Single Lumen Midline Seldinger Kit	5	20

#### 4FR Kit Range

Catalogue No.	Description	Box Qty	Carton Qty
A130S-04106	4Fr x 06cm Single Lumen Midline Seldinger Kit	5	20
A130S-04108	4Fr x 08cm Single Lumen Midline Seldinger Kit	5	20
A130S-04110	4Fr x 10cm Single Lumen Midline Seldinger Kit	5	20



## Standard Kit Components\*

- Catheter
- Safety Introducer Needle
- Guidewire
- Securement Device
- Wristband

\*Components pictured may vary based on availability

Also available in Maximum Barrier kit configurations.

For ordering information or for questions about the Health Line ProV Midline ST, please contact your Health Line Sales Consultant at (877) 847-4542 or email customerservice@hlic.net



