

SYNERGY-XS CT PICC™


Smaller Vein PICC Access
with Power Injectability



HEALTH LINE
MEDICAL PRODUCTS

Health Line now offers a complete line of 3Fr Single and 4Fr Dual Lumen power injectable Synergy CT PICC™ catheters for patients with smaller veins. These XS CT PICCs offer the same ease and effectiveness as our popular Synergy CT PICC™ line with the advantage of giving clinicians more options in choosing the best PICC size for the patient.

The Synergy-XS CT PICCs™ are available in a variety of kit configurations including basic, nurses, and Interventional Radiology to give clinicians an assortment of options to meet patient needs.



Features and Benefits

One Piece Cannula Construction - No bonds or welds helps to ensure catheter will not break.

CVP Monitoring - SYNERGY CT PICC™ (power PICCs) are indicated for central vein pressure monitoring.

Polyurethane Catheter Construction - Provides radiopacity, durability, and strength.


Low Profile Hub - Allows for better patient comfort and easy securement.

Grip-Lok™ - Stabilizes and locks catheter securely in place reducing needlestick injuries.

Safety Components - Kit trays include a range of safety components such as, safety scalpel and safety introducer needle to improve patient and clinician safety.

Synergy-XS CT PICC™ (Peripherally Inserted Central Catheters)

Catalogue No.	PICC Description	Guide Wire	Introducer	Box Qty.	Carton Qty.
NURSES KIT					
A14-03160-N	3f x 60cm Single Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	3.5f x 05cm PTFE Tearaway	5	20
A14-04260-N	4f x 60cm Dual Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	4.5f x 05cm PTFE Tearaway	5	20
BASIC KIT					
A14-03160	3f x 60cm Single Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	3.5f x 10cm PTFE Tearaway	5	20
A14-04260	4f x 60cm Dual Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
IR KIT					
A14-03160-IR	3f x 60cm Single Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	3.5f x 10cm PTFE Tearaway	5	20
A14-04260-IR	4f x 60cm Dual Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20



Standard Kit Components

Component	Quantity	Component	Quantity
Synergy-XS CT PICC™	1	Safety Scalpel #11	1
Needle Free Valve	1 or 2	24" Paper Tape Measure	1
Twisted Wire Stylet and Adaptor w/ Sideport	1	PICC Stylet TAG	1
21GA X 7CM Safety Seld.Echo.Needle	1	HLIC PICC Patient Info Guide	1
PTFE Tearaway Introducer	1	HLIC CT PICC Chart Sticker	1
Nitinol Straight Mandrel GW w/ Tungsten Tip	1	HLIC CT PICC IFU	1
10cc Disposable Luer Lock Syr	1	Securment Device	1

For ordering information or for questions about the Synergy-XS CT PICC™, please contact your Health Line Sales Consultant at (877) 847-4542 and visit us at www.hlic.net.





Health Line International Corp.
5675 West 300 South
Salt Lake City, Utah 84104

Phone: (877) 847-4542
Fax: (855) 228-1336
www.hlic.net

ISO13485
Certified Quality System
MKT-100-08 Rev. 02