


# SYNERGY™

## PICC Access with Power Injectability



The SYNERGY™ CT PICC offers the ease and effectiveness of PICC access with power injectability, creating a SYNERGY of intravenous therapy options for healthcare practitioners.



## Features and Benefits

**Maximum Injection Rates** – Allows for the injection of contrast media at rates up to 5ml/sec for CECT Scans.

**CVP Monitoring** – SYNERGY™ power PICCs are indicated for central vein pressure monitoring.

**Polyurethane Construction** – Reverse taper technology helps reduce post insertion bleeding while maintaining radiopacity and patient comfort.

**Standardized Trays and Components** – Trays are organized with familiar industry components for clinician safety and user friendly procedural steps; including tearaway sheaths, echogenic safety needles, and premium guide wires.

**Low Profile Hub** – Allows for better patient comfort and easy securement.

**Multiple Kit Configurations** – Available in Nursing, Basic, and I.R. Long Wire tray configurations.

## Synergy™ CT Peripherally Inserted Central Catheters (PICC)

	Catalogue No.	PICC Description	Guide Wire	Introducer	Box Qty.	Carton Qty.
NURSES KIT	<b>A14-04160-N</b>	4f x 60cm Single Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05160-N</b>	5f x 60cm Single Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05260-N</b>	5f x 60cm Dual Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06260-N</b>	6f x 60cm Dual Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06360-N</b>	6f x 60cm Triple Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
BASIC KIT	<b>A14-04160</b>	4f x 60cm Single Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05160</b>	5f x 60cm Single Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05260</b>	5f x 60cm Dual Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06260</b>	6f x 60cm Dual Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06360</b>	6f x 60cm Triple Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
IR KIT	<b>A14-04160-IR</b>	4f x 60cm Single Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05160-IR</b>	5f x 60cm Single Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05260-IR</b>	5f x 60cm Dual Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06260-IR</b>	6f x 60cm Dual Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06360-IR</b>	6f x 60cm Triple Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20

## Standard Tray Components

SYNERGY™ CT PICC, Needle Free Valve (1-3), Pre-Loaded Twisted Wire Stylet with Flushable T-Port, Safety Seldinger Echo Tip Needle, PTFE Tearaway Introducer, Nitinol Straight Mandrel Guide Wire, 10cc Disposable Luer Lock Syringe, #11 Safety Scalpel, 24" Paper Tape Measure, PICC Stylet TAG, PICC Patient Info Guide, PICC Chart Sticker, IFU, Securement Device.

**For ordering information or for questions about the SYNERGY™ PICC Line, please contact your Health Line Sales Consultant at (877) 847-4542.**



Health Line International Corp.  
5675 West 300 South  
Salt Lake City, Utah 84104

Phone: (877) 847-4542  
Fax: (855) 228-1336  
[www.hlic.net](http://www.hlic.net)

**ISO13485**  
**Certified Quality System**  
**MKT-300-02 Rev. 01**