

# SYNERGY CT PICC™


PICC Access with  
Power Injectability



HEALTH LINE  
MEDICAL PRODUCTS

PART OF  NUMANTEC®  
HUMAN INNOVATION TECHNOLOGY

The SYNERGY CT PICC™ offers the ease and effectiveness of PICC access with power injectability, creating a SYNERGY of intravenous therapy options for healthcare practitioners.



## Features and Benefits

**Maximum Injection Rates** - Allows for the injection of contrast media at rates up to 5ml/sec for CECT Scans.

**CVP Monitoring** - SYNERGY CT PICC™ (power PICCs) are indicated for central vein pressure monitoring.

**Polyurethane Construction** - Reverse taper technology helps reduce post insertion bleeding while maintaining radiopacity and patient comfort.

**Standardized Trays and Components** - Trays are organized with familiar industry components for clinician safety and user friendly procedural steps; including tearaway sheaths, echogenic safety needles, and premium guide wires.

**Low Profile Hub** - Allows for better patient comfort and easy securement.

**Multiple Kit Configurations** - Available in Nursing, Basic, and I.R. Long Wire tray configurations.

## Synergy CT PICC™ (Peripherally Inserted Central Catheters)

	Catalogue No.	PICC Description	Guide Wire	Introducer	Box Qty.	Carton Qty.
<b>NURSES KIT</b>	<b>A14-04160-N</b>	4f x 60cm Single Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	4.5f x 05cm PTFE Tearaway	5	20
	<b>A14-05160-N</b>	5f x 60cm Single Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	5.5f x 05cm PTFE Tearaway	5	20
	<b>A14-05260-N</b>	5f x 60cm Dual Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	5.5f x 05cm PTFE Tearaway	5	20
	<b>A14-06260-N</b>	6f x 60cm Dual Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	6.5f x 05cm PTFE Tearaway	5	20
	<b>A14-06360-N</b>	6f x 60cm Triple Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	6.5f x 05cm PTFE Tearaway	5	20
<b>BASIC KIT</b>	<b>A14-04160</b>	4f x 60cm Single Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05160</b>	5f x 60cm Single Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05260</b>	5f x 60cm Dual Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06260</b>	6f x 60cm Dual Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06360</b>	6f x 60cm Triple Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
<b>IR KIT</b>	<b>A14-04160-IR</b>	4f x 60cm Single Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05160-IR</b>	5f x 60cm Single Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-05260-IR</b>	5f x 60cm Dual Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06260-IR</b>	6f x 60cm Dual Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	<b>A14-06360-IR</b>	6f x 60cm Triple Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20

### Standard Tray Components

SYNERGY CT PICC™, Needle Free Valve (1-3), Pre-Loaded Twisted Wire Stylet with Flushable T-Port, Safety Seldinger Echo Tip Needle, PTFE Tearaway Introducer, Nitinol Straight Mandrel Guide Wire, 10cc Disposable Luer Lock Syringe, #11 Safety Scalpel, 24" Paper Tape Measure, PICC Stylet TAG, PICC Patient Info Guide, PICC Chart Sticker, IFU, Securement Device.

**For ordering information or for questions about the SYNERGY CT PICC™ Line, please contact your Health Line Sales Consultant at (877) 847-4542.**



