


SYNERGY CT PICC™

PICC Access with Power
Injectability



The SYNERGY CT PICC™ offers the ease and effectiveness of PICC access with power injectability, creating a SYNERGY of intravenous therapy options for healthcare practitioners.



Features and Benefits

Maximum Injection Rates: Allows for the injection of contrast media at rates up to 5ml/sec for CECT Scans.

CVP Monitoring: SYNERGY CT PICC™ (power PICCs) are indicated for central vein pressure monitoring.

Polyurethane Construction: Reverse taper technology helps reduce post insertion bleeding while maintaining radiopacity and patient comfort.

Standardized Trays and Components: Trays are organized with familiar industry components for clinician safety and user friendly procedural steps; including tearaway sheaths, echogenic safety needles, and premium guide wires.

Low Profile Hub: Allows for better patient comfort and easy securement.

Multiple Kit Configurations: Available in Nursing, Basic, and I.R. Long Wire tray configurations

Synergy CT PICC™ (Peripherally Inserted Central Catheters)

	Catalogue No.	PICC Description	Guide Wire	Introducer	Box Qty.	Carton Qty.
NURSES KIT	A14-04160-N	4f x 60cm Single Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	4.5f x 05cm PTFE Tearaway	5	20
	A14-05160-N	5f x 60cm Single Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	5.5f x 05cm PTFE Tearaway	5	20
	A14-05260-N	5f x 60cm Dual Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	5.5f x 05cm PTFE Tearaway	5	20
	A14-06260-N	6f x 60cm Dual Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	6.5f x 05cm PTFE Tearaway	5	20
	A14-05360-N	5f x 60cm Triple Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	6.5f x 05cm PTFE Tearaway	5	20
	A14-06360-N	6f x 60cm Triple Lumen CT PICC	45cm Nitinol GW W/ Tungsten Tip	6.5f x 05cm PTFE Tearaway	5	20
BASIC KIT	A14-04160	4f x 60cm Single Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
	A14-05160	5f x 60cm Single Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	A14-05260	5f x 60cm Dual Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	A14-06260	6f x 60cm Dual Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	A14-05360	5f x 60cm Triple Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	A14-06360	6f x 60cm Triple Lumen CT PICC	70cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
IR KIT	A14-04160-IR	4f x 60cm Single Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	4.5f x 10cm PTFE Tearaway	5	20
	A14-05160-IR	5f x 60cm Single Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	A14-05260-IR	5f x 60cm Dual Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	5.5f x 10cm PTFE Tearaway	5	20
	A14-06260-IR	6f x 60cm Dual Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	A14-05360-IR	5f x 60cm Triple Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20
	A14-06360-IR	6f x 60cm Triple Lumen CT PICC	130cm Nitinol GW W/ Tungsten Tip	6.5f x 10cm PTFE Tearaway	5	20

Standard Tray Components

SYNERGY CT PICC™, Needle Free Valve (1-3), Pre-Loaded Twisted Wire Stylet with Flushable T-Port, Safety Seldinger Echo Tip Needle, PTFE Tearaway Introducer, Nitinol Straight Mandrel Guide Wire, 10cc Disposable Luer Lock Syringe, #11 Safety Scalpel, 24" Paper Tape Measure, PICC Stylet TAG, PICC Patient Info Guide, PICC Chart Sticker, IFU, Securement Device.

For ordering information or for questions about the SYNERGY CT PICC™ Line, please contact your Health Line Sales Consultant at (877) 847-4542.





HEALTH LINE
MEDICAL PRODUCTS

PART OF



NUMANTEC®
HUMAN INNOVATION TECHNOLOGY

Health Line International Corp.
260 North Ace Yeager Court, Unit D
Salt Lake City, Utah 84116

Phone: (877) 847-4542
customerservice@hlic.net
www.hlic.net

ISO13485
Certified Quality System
MKT-300-02 Rev. 06